

DMT-301 Tan Dental Modeling

General Information

DMT-301 is formulated to work with DLP and SLA 3D printers using 385-405 nm light. It is specifically tuned to print in open mode on the Form2 3D printer. DMT-301 will extend the window lifetime of the windows of printers that use PDMS interfaces. This resin reliably produces prints with low shrinkage as well as good positive and negative resolution.

Chemical Data

		Method
Viscosity (cps @25 °C)	710	Brookfield SP #31
	Post-cured	
Tensile Modulus (mPa)	1270	ASTM D 638-1V
Tensile Strength (mPa)	58	ASTM D 638-1V
Elongation (%)	5	ASTM D 638-1V
Shore Hardness	86.5 D	

**Parts were post-cured for 10min a side with a 395-405nm LED at 10mW/cm2.*

**Properties may vary in orientation and post-treatment.*

Form2 Settings

DMT-301 is recommended to print with Grey V4 settings in Open Mode.

Post Cure Procedure

Parts should be post cured under a either a 395-405 nm LED for 10 minutes on a side at a minimum of 10 mW/cm² or a broadband UV light source at 10 minutes a side to bring the parts up to their full properties

Storage and Handling

DMT-301 stable for up to 1 year at ambient conditions.

This product is light sensitive and should not be exposed to daylight, UV light and artificial lighting during storage.

Exposure to daylight, UV light and artificial lighting should be kept to a minimum during handling.

Arkema Inc., 3D Printing Development Center

1880 S. Flatiron Ct STE J, Boulder, CO 80301

Please consult Arkema's disclaimer regarding the use of Arkema's products on <https://www.arkema.com/global/en/products/product-safety/disclaimer/>

Arkema - a French société anonyme, registered in the Nanterre (France) Trade and Companies Register under the number 319 632 790

[arkema.com](https://www.arkema.com)

ARKEMA